

# Stereo Zoom Microscope 10-40x

---

Instruction Manual for  
MI-40STER – Stereo Zoom Microscope 10-40x

---

The logo for Home Science Tools features the words "HOME SCIENCE TOOLS" in a bold, sans-serif font. Above the letter "I" in "SCIENCE" are two small circles connected by a horizontal line, resembling a pair of glasses or a microscope. Below the main text, the tagline "THE GATEWAY TO DISCOVERY" is written in a smaller, all-caps, sans-serif font.

**HOME SCIENCE TOOLS**  
THE GATEWAY TO DISCOVERY

1.800.860.6272  
[www.homesciencetools.com](http://www.homesciencetools.com)

Copyright 2005 by Home Training Tools, Ltd.  
All rights reserved.

Welcome to an exciting world of discovery with your new **Stereo Zoom Microscope!** This low-power microscope is designed for viewing whole objects, such as flowers, rocks, or insects. This manual will give you a familiarity with the different features of your microscope, how to use them, and how to preserve your investment by proper maintenance and care.


motion. Repeat with a second paper moistened with solution if necessary. Repeat once again with a piece of dry lens paper until the lens is clean and dry. **Do not spray lens cleaner directly on the lens.**

## Table of Contents

General Microscope Care.....	2
Unpacking .....	2
Cleaning .....	2
Features & Definitions .....	2
Microscope Diagram .....	2
Description of Components.....	2
Operating Procedure .....	3
Getting Started.....	3
Using the Binocular Head .....	4
Maintenance .....	4
Changing the Top Bulb .....	4
Changing the Bottom Bulb .....	4
Adjusting Tension.....	4
Specifications .....	5
Troubleshooting .....	5
Warranty .....	5
Ideas for Using Your Microscope .....	6
Microscope Observation Worksheet.....	7

## Features & Definitions

### Microscope Diagram



## General Microscope Care

### Unpacking

Your stereo microscope is shipped in a two-part Styrofoam case. Keep this case for storage, transport, and shipping. It is perfect packing material should you ever need to send your microscope in for repairs covered by the warranty.


Avoid touching the lens surfaces on the eyepiece or objective lens, as finger prints will decrease image quality.

### Cleaning

The best optical quality can be compromised by dirty lenses. Using a dustcover and cleaning the lenses regularly will greatly enhance your microscope use.

To clean lens surfaces, remove dust by using a soft brush or a can of compressed air. Then moisten a piece of lens paper (our item MI-PAPER) with some lens cleaning solution (MI-LENSCLN). Gently clean the eyepiece and objective lens exterior surface using a circular

### Description of Components

- 1. Binocular eyepieces:** This is the part of the microscope that you look through. The binocular eyepieces contain lenses that magnify 10x and provide an unreversed 3D image. They are inclined at an angle for comfortable viewing.
- 2. Stereo head:** The head rotates 360° so that multiple users can look in the eyepieces comfortably without moving the microscope itself (to rotate, loosen the head set screw on the left side of the microscope). 
- 3. Diopter:** This knurled band on each eyepiece is used to adjust the focus for differences between your eyes. Instructions for doing this are on page 4.

4. **Zoom knob:** This knob is used to change magnification. It allows you to “zoom” from 10x to 40x.
5. **Objective turret:** This turret contains the lenses closest to the specimen. The objective lenses have magnification between 1x and 4x (providing a total magnification of 10x-40x when multiplied with the 10x of the eyepieces).
6. **Locking knob:** The binocular head is mounted on a post and can be raised, lowered, or turned around by loosening the locking knob on the back of the post.
7. **Top illuminator:** This bulb-holder holds the 10-watt halogen bulb that shines down on the specimen. Use this light when your specimen is opaque or solid (when light cannot pass through it from below).
8. **Focus knob:** This knob is used to raise or lower the objective lens until the image is in focus.
9. **Stage:** The stage is the platform that supports the specimen below the objective lens.
10. **Head stop:** This sets the lowest position the head can drop. For normal use it can be left in the lowest position. If you are examining tall specimens, adjust the ring so that the head cannot hit the specimen.
11. **Stage plate screw:** The stage plate can be removed and changed by loosening this screw.
12. **Stage clips:** These clips can be used to hold thin specimens in place.
13. **Stage plate:** This microscope comes with two stage plates. The glass plate is used with bottom lighting, and the reversible black/white plate is used with top lighting to help you get the best contrast.
14. **Bottom illuminator:** Another 10-watt halogen bulb is located beneath the stage plate. Use this light for translucent specimens.
15. **Illuminator control:** This allows you to choose three different light settings: top lighting, bottom lighting, or top and bottom together.

## Operating Procedure

Now that you have an overview of what each component of your microscope is for, you can

follow this step-by-step procedure to help you get started using it.

## Getting Started

1. Set your microscope on a tabletop or other flat sturdy surface where you will have plenty of room to work. Plug the microscope's power cord into an outlet, making sure that the excess cord is out of the way so no one can trip over it or pull it off of the table.
2. Flip the switch to turn on your microscope's light source. Use top lighting for opaque specimens and bottom lighting for translucent specimens. Some specimens have both opaque and translucent parts. For these use top and bottom lighting together. **Warning:** The top light can get very hot. Use care touching the top light housing during use.
3. Center a specimen on the stage plate. If you are using top lighting, insert the reversible black/white stage plate (use the dark side if the specimen is light colored). To change or reverse the plate, loosen the stage plate screw until you are able to pop the plate out. Turn the plate over and tighten the screw to lock it in place.
4. If your specimen is thin and flat, or if its edges curl up easily, use the stage clips to hold it in place. To do this, pull up the pointed end of one stage clip and slide it over one end of the specimen, then do the same with the stage clip on the other side. If your specimen is larger than the stage plate, turn the stage clips out so that they hang off the stage; this will give you more room to work.
5. You may need to adjust the height of the head in order to find a good working distance between the specimen and the objective lens. Do this by loosening the locking knob, moving the head to the appropriate position, and tightening the locking knob.
6. Turn the zoom knob away from you until the microscope is on its lowest power.
7. Slowly turn the focus knob until the specimen comes into view. Once you can see the outline of the specimen, turn the knob even more slowly until it is focused as sharply as possible. Once you have focused your specimen, you can move it around to see other parts of it. You may need to refocus slightly on each new area. **Note:** with this microscope you will often be viewing three-dimensional specimens that have many

different levels. You will not be able to focus every feature clearly at the same time.

8. For higher magnification, turn the zoom knob toward you. The image will move slightly out of focus as you zoom, so you will need to make minor focus adjustments each time you stop.

### Using the Binocular Head

To use the binocular head to the most advantage, you must set the interpupillary distance to match the distance between the pupils of your eyes. You must also adjust the diopter to compensate for focusing differences between your two eyes. Each user of the microscope must make these adjustments for his or her own eyes. To do so, follow these steps:

1. Focus the microscope on a small specimen in the center of the stage plate.
2. Focus your eyes on the specimen.
3. Pull your eyes back from the eyepieces about 1". In your peripheral vision you will see two field view circles overlapping each other.
4. Open or close the distance between the eyepieces by twisting them apart or pushing them together until the two circles merge together and appear as one circle. The interpupillary distance is set correctly when you see just one field view circle.
5. Close your left eye and turn the knurled diopter band until the specimen is in focus for your right eye.
6. Close your right eye and bring the specimen into sharp focus for your left eye by turning the knurled diopter band on the left eyepiece.

## Maintenance

### Changing the Top Bulb

If your top microscope bulb burns out, follow these steps to replace it:

1. Obtain the correct 10-watt halogen replacement bulb (our item MI-BULB9) with reflector.
2. Unplug your microscope from the power supply and allow it to cool before replacing the bulb.
3. Unscrew the cap that covers the bulb.
4. Grasp the bulb and pull it straight out from the socket.



Note: The glass reflector is part of the bulb and comes out with it.

5. Use a tissue or cloth to grasp the new bulb and insert it into the socket.

### Changing the Bottom Bulb


If your bottom microscope bulb burns out, follow these steps to replace it:

1. Obtain the correct 10-watt halogen replacement bulb (our item MI-BULB8).
2. Unplug your microscope from the power supply and allow it to cool before replacing the bulb.
3. Remove the stage plate by loosening the stage plate screw until the plate can be removed easily. Remove the blue filter located under the stage plate.
4. Pull the bulb straight out of its socket.
5. Use a tissue or cloth to grasp the new bulb and insert it into the socket.
6. Replace the stage plate and secure it in place with the stage plate screw.



### Adjusting Tension

The focus tension is pre-adjusted by the manufacturer, but if it falls out of adjustment, the microscope head will drift down under its own weight and the image will move out of focus. The tension adjustment collar is located between the microscope arm and the focus knob on the left side (when the stage is facing you). To adjust the tension, follow these steps:



1. Turn the collar clockwise to tighten or counter-clockwise to loosen. (Use a rubber band to help you grasp the collar, or a small C-wrench.)
2. Tighten only enough to keep the stage from drifting downward.

## Specifications

- Premium widefield 10x eyepieces, fully coated
- Inclined 45° head rotates 360° with 55-75mm interpupillary distance and one diopter (+/- 5 diopters)
- Fully coated 1x-4x zoom objective provides magnification from 10x to 40x
- 80mm working distance, extra-large 23mm field-of-view at 10x, 6.5 mm field-of-view at 40x
- All metal rack-and-pinion focusing, focus knob has tension adjustment
- Stage with stage clips, frosted glass stage plate, and reversible white/black stage plate
- Maximum specimen height of 110mm on-stage
- Halogen 10-watt, 12 volt top and bottom light illumination with grounded 110 volt cord
- Operates with top only, bottom only, or top and bottom illumination together

## Troubleshooting

If you are experiencing difficulty with your microscope, try these troubleshooting techniques:

Problem	Possible Reason and Solution
<b>Light fails to operate</b>	<ol style="list-style-type: none"><li>1. The AC power cord is not connected. <i>Connect the cord to an outlet.</i></li><li>2. The bulb is burned out. <i>Replace the bulb. (See "Changing the Bulb," p. 4.)</i></li><li>3. An incorrect bulb is installed. <i>Replace with the correct bulb.</i></li></ol>
<b>Light flickers</b>	<ol style="list-style-type: none"><li>1. The bulb is not properly inserted into the socket. <i>Properly insert the bulb.</i></li><li>2. The bulb is about to burn out. <i>Replace the bulb.</i></li></ol>
<b>Poor resolution, image not sharp, spots in the field</b>	<ol style="list-style-type: none"><li>1. The objective or eyepiece lenses are dirty. <i>Clean the lenses. (See "Cleaning," p. 2.)</i></li></ol>

## Warranty

Home Science Tools warrants this microscope to be free from defects in material and workmanship under normal use and service for the life of the instrument. This warranty does not cover light bulbs, batteries, or damage due to misuse, abuse, alterations, or accident. Warranty does not cover lenses that have become inoperable due to excessive dirtiness as a result of misuse or lack of normal maintenance.

Any cameras and software supplied with this microscope are warranted for one year from the date of purchase.

You will need to return your microscope freight prepaid for warranty service to Home Science Tools, or the repair facility we designate. We will repair or replace your microscope at no charge and return it freight prepaid to you. Please call 1-800-860-6272 to arrange warranty service before returning this instrument.

Please note that warranties apply only to the original purchaser and are not transferable.

## Ideas for Using Your Microscope

Your stereo microscope is a versatile instrument than can be used to view a variety of specimens. This section contains various suggestions for what to study.

Clear plastic or glass petri dishes are great for viewing live or messy objects with a stereo microscope because they fit well on the stage and keep everything adequately contained. The suggestions below are just a few things you can view with petri dishes. Place the item or items to be viewed in the bottom of a petri dish and position it on the stage plate of your microscope. Use top or bottom lighting.

### **Observe the habits of live insects.**

Collect insects in the bottom of a petri dish and cover with its lid to keep insects from escaping. Be careful not to leave the light source shining on the insects for too long as the heat could eventually kill them.

### **Study a shallow dish of pond water, daphnia, or fairy shrimp.**

Watch them closely as these tiny creatures swim, dive, and eat.

### **Examine a soil sample to see the different materials that comprise it.**

Soils with a lot of sand or clay are particularly interesting to look at. You might even want to collect soil samples from several different spots and compare and contrast what you see in each sample.

### **Dissect a flower to learn about the beauty and intricacies of all its parts.**

Carefully pull the flower petals and inside parts off of the stem trying not to damage or tear them. See if you can identify the parts using a flower identification book. Stick one or two of the parts on your microscope to get a closer look. If there was a lot of pollen on the flower, try putting the pollinated parts, or loose pollen, into a petri dish and check it out with your microscope. *(Note: This is not a good experiment to do if you have bad allergies!)*

### **Compare the types of minerals and crystals in different rock specimens.**

You can break off small pieces of larger rocks by knocking them together or using a rock pick. Put any small shards or pieces of the broken rocks into a petri dish for easy viewing.

### **Make a simple prepared slide.**

To make a slide, tear a 2½-3" long piece of Scotch tape and set it sticky side up on the kitchen table or other work area. Fold over about ½" of the tape on each end to form finger holds on the sides of the slide. Next, sprinkle a few grains of salt, sugar, ground coffee, or sand in the middle of the sticky part of the slide. Carefully observe the differences between different grains.

Hair and thread also work well on homemade tape slides. Collect samples of hair from family members or pets and stick one hair from each sample on a tape slide. Label each slide and view them one at a time with your microscope. Write down your observations about each to see how hairs from humans and animals differ. You can also look at threads or fibers from furniture, rugs or clothing from around your house.

### **Record your observations.**

In the field of science, recording observations while performing an experiment is one of the most useful tools available. Early scientists often kept very detailed journals of the experiments they performed, making entries for each individual experiment and writing down virtually everything they saw. These entries often included drawings and detailed descriptions as well as the procedures they used, the data they collected, and conclusions drawn from their experimentation.

Our Microscope Observation worksheet (on the next page) will help you keep track of the things that you study with your microscope and remember what you have learned. Blanks are provided for recording general information about each specimen, such as its type and the date it was collected. In addition, there is space to write down your observations and make sketches of what you see.

# Microscope Observation Worksheet

Name of specimen: \_\_\_\_\_

Date specimen was collected: \_\_\_\_\_

Collected from: \_\_\_\_\_

## Observations

---


---

---


---

---


## Sketches




Lowest power



Highest power



Other: \_\_\_\_\_



Other: \_\_\_\_\_